


GAD67 Rabbit mAb

ABZ4852

Product Overview

Name	GAD67 Rabbit mAb
Accession(Primary)	Q99259
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 67 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q99259)

Gene Name GAD1**About Gene Symbol**

Western blot analysis of GAD67 in (1) Jurkat lysates; (2) C6 lysates using GAD67 antibody.

