


SynGAP Rabbit mAb

ABZ4854

Product Overview

Name	SynGAP Rabbit mAb
Accession(Primary)	Q96PV0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 148 kDa; Observed MW: 148 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q96PV0)

Gene Name SYNGAP1**About Gene Symbol**

Western blot analysis of SynGAP in Human fetal brain lysates using SynGAP antibody.

