

ATG4C Rabbit mAb


ABZ4856

Product Overview

Name	ATG4C Rabbit mAb
Accession(Primary)	Q96DT6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:50-1:100

Target information (Q96DT6)

Gene Name	ATG4C
About Gene Symbol	ATG4C (Autophagy Related 4C Cysteine Peptidase) is a protease enzyme involved in autophagy. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify



Western blot analysis of Atg4C in Jurkat lysates
using ATG4C antibody.