


## ALIX Rabbit mAb

ABZ4860

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | ALIX Rabbit mAb                            |
| <b>Accession(Primary)</b>   | Q8WUM4                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 96 kDa; Observed MW: 96 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                           |

### Target information (Q8WUM4)

**Gene Name** PDCD6IP**About Gene Symbol**

Western blot analysis of ALIX in K562 lysates using ALIX antibody.

