

Podoplanin Rabbit mAb

ABZ4862


Product Overview

Name	Podoplanin Rabbit mAb
Accession(Primary)	Q86YL7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q86YL7)

Gene Name PDPN

About Gene Symbol



Western blot analysis of Podoplanin in Human placenta lysates using Podoplanin antibody.

