

# Podoplanin Rabbit mAb

ABZ4863


## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Podoplanin Rabbit mAb                      |
| <b>Accession(Primary)</b>   | Q86YL7                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 17 kDa; Observed MW: 36 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500        |

## Target information (Q86YL7)

**Gene Name** PDPN

**About Gene Symbol**



Western blot analysis of Podoplanin in 293T lysates using Podoplanin antibody.

