

CNOT7 Rabbit mAb


ABZ4866

Product Overview


Name	CNOT7 Rabbit mAb
Accession(Primary)	Q60809
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:500-1:1000; FC-1:200-1:1000; IP-1:20-1:50; IHC-1:50-1:100

Target information (Q60809)


Gene Name	Cnot7
About Gene Symbol	



Western blot analysis of CNOT7 in F9, C6, PC-12, 3T3, Raw264.7 lysates using CNOT7 antibody.



Western blot analysis of CNOT7 in Molt4, Raji, COS7, HeLa lysates using CNOT7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using CNOT7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.