

GIV Rabbit mAb


ABZ4867

Product Overview


Name	GIV Rabbit mAb
Accession(Primary)	Q3V6T2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 216 kDa; Observed MW: 216 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:500-1:2000; FC-1:50-1:100; IP-1:20-1:50

Target information (Q3V6T2)

Gene Name	CCDC88A
About Gene Symbol	



Western blot analysis of GIV in mouse brain, Raw264.7 lysates using GIV antibody.



Western blot analysis of Girdin in HeLa, A549, HL-60, U2OS, C6 lysates using Girdin antibody.