

Fascin Rabbit mAb



ABZ4869

Product Overview

Name	Fascin Rabbit mAb
Accession(Primary)	Q16658
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:100-1:1000; IP-1:10-1:100


Target information (Q16658)

Gene Name	FSCN1
About Gene Symbol	




Immunocytochemistry analysis of Fascin (green) in 3T3 using Fascin antibody, and DAPI (blue).

Western blot analysis of Fascin in mouse skin, rat skin lysates using Fascin antibody.



Western blot analysis of Fascin in HeLa, A549, HepG2, COS7 lysates using Fascin antibody



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Fascin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.