

NRF1 Rabbit mAb


ABZ4870

Product Overview


Name	NRF1 Rabbit mAb
Accession(Primary)	Q16656
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; FC-1:10-1:100; IP-1:10-1:100

Target information (Q16656)

Gene Name	NRF1
About Gene Symbol	



Western blot analysis of Nrf1 in HeLa, U87-MG and Jurkat lysates using Nrf1 antibody.



Western blot analysis of Nrf1 in K562, Ramos and Jurkat lysates using NRF1 antibody.