

NEDD8 Rabbit mAb


ABZ4874

Product Overview

Name	NEDD8 Rabbit mAb
Accession(Primary)	Q15843
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 9 kDa; Observed MW: 10,90 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:100-1:200; FC-1:200-1:500; IP-1:20-1:50

Target information (Q15843)

Gene Name	NEDD8
About Gene Symbol	



Western blot analysis of NEDD8 in K562 lysates
using NEDD8 antibody.