

Smad2/3 Rabbit mAb


ABZ4875

Product Overview


Name	Smad2/3 Rabbit mAb
Accession(Primary)	Q15796
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52,60 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:100-1:200; FC-1:200-1:500; IP-1:20-1:50

Target information (Q15796)

Gene Name	SMAD2
About Gene Symbol	SMAD2 (SMAD Family Member 2) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Smad2/3 (green) in HT-1080 using Smad2/3 antibody, and DAPI (blue).



Western blot analysis of Smad2/3 in HeLa, 293, C6, C2C12 lysates using Smad2/3 antibody.