

IKB beta Rabbit mAb


ABZ4877

Product Overview


Name	IKB beta Rabbit mAb
Accession(Primary)	Q15653
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:20-1:50; FC-1:20-1:50; IP-1:10-1:100

Target information (Q15653)

Gene Name	NFKB1B
About Gene Symbol	



Western blot analysis of IKB beta in Jurkat lysates using IKB beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human uterus using IKB beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.