

TOMM20 Rabbit mAb


ABZ4878

Product Overview


Name	TOMM20 Rabbit mAb
Accession(Primary)	Q15388
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (Q15388)


Gene Name	TOMM20
About Gene Symbol	



Western blot analysis of TOMM20 in HepG2 lysates using TOMM20 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using TOMM20 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of TOMM20 in HeLa using TOMM20 antibody.

