

TOMM20 Rabbit mAb


ABZ4879

Product Overview


Name	TOMM20 Rabbit mAb
Accession(Primary)	Q15388
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (Q15388)


Gene Name	TOMM20
About Gene Symbol	




Western blot analysis of TOMM20 in 3T3 lysates using TOMM20 antibody.



Western blot analysis of TOMM20 in MCF-7, T47D lysates using TOMM20 antibody.



Immunocytochemistry analysis of TOMM20(green)
in HeLa using TOMM20 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-
embedded Human breast cancer using TOMM20
antibody. High-pressure and temperature Sodium
Citrate pH 6.0 was used for antigen retrieval.