


PKC mu/PKD Rabbit mAb

ABZ4880

Product Overview

Name	PKC mu/PKD Rabbit mAb
Accession(Primary)	Q15139
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q15139)

Gene Name PRKD1**About Gene Symbol**

Western blot analysis of PKC mu in HeLa lysates using PKD1 antibody.

