

Lipin 1 Rabbit mAb


ABZ4881

Product Overview

Name	Lipin 1 Rabbit mAb
Accession(Primary)	Q14693
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 99 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:10-1:100; IP-1:20

Target information (Q14693)

Gene Name	LPIN1
About Gene Symbol	



Western blot analysis of Lipin 1 in Jurkat lysates using Lipin 1 antibody.