

TRIM25 Rabbit mAb


ABZ4882

Product Overview

Name	TRIM25 Rabbit mAb
Accession(Primary)	Q14258
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 71 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:1000; FC-1:50-1:100; IP-1:20-1:100

Target information (Q14258)

Gene Name	TRIM25
About Gene Symbol	TRIM25 (Tripartite Motif Containing 25) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and i



Western blot analysis of TRIM25 in HeLa lysates
using TRIM25 antibody.