

## IKK epsilon Rabbit mAb


ABZ4883

### Product Overview

|                             |                                                         |
|-----------------------------|---------------------------------------------------------|
| <b>Name</b>                 | IKK epsilon Rabbit mAb                                  |
| <b>Accession(Primary)</b>   | Q14164                                                  |
| <b>Host Species</b>         | Rabbit                                                  |
| <b>Clonality</b>            | Monoclonal                                              |
| <b>Isotype</b>              | IgG                                                     |
| <b>Species Reactivity</b>   | Human                                                   |
| <b>Calculated MW</b>        | Calculated MW: 80 kDa; Observed MW: 80 kDa              |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>FC-1:200-1:1000 |

### Target information (Q14164)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | IKBKE |
| <b>About Gene Symbol</b> |       |



Western blot analysis of IKKi/IKKe in HeLa, Raji, Ramos, MCF-7 lysates using IKK epsilon antibody.