

Cullin 4A Rabbit mAb


ABZ4884

Product Overview

Name	Cullin 4A Rabbit mAb
Accession(Primary)	Q13619
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:20-1:50

Target information (Q13619)

Gene Name	CUL4A
About Gene Symbol	CUL4A (Cullin 4A) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cullin 4a in HeLa lysates
using Cullin 4A antibody.