


Cullin 3 Rabbit mAb

ABZ4885

Product Overview

Name	Cullin 3 Rabbit mAb
Accession(Primary)	Q13618
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 89 kDa; Observed MW: 89 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:100

Target information (Q13618)

Gene Name CUL3**About Gene Symbol**

Western blot analysis of Cullin 3 in SH-SY5Y
lysates using Cullin 3 antibody.

