

Stromal Interaction Molecule 1 Rabbit mAb


ABZ4886

Product Overview

Name	Stromal Interaction Molecule 1 Rabbit mAb
Accession(Primary)	Q13586
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (Q13586)

Gene Name	STIM1
About Gene Symbol	STIM1 (Stromal Interaction Molecule 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Stromal interaction molecule 1 in HeLa lysates using Stromal Interaction Molecule 1 antibody.