

CaMKII Rabbit mAb


ABZ4887

Product Overview

Name	CaMKII Rabbit mAb
Accession(Primary)	Q13554
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC:1:200-1:500; ICC/IF-1:20-1:100; FC-1:20-1:100

Target information (Q13554)

Gene Name	CAMK2B
About Gene Symbol	



Western blot analysis of CaMKII in HeLa lysates
using CaMKII beta antibody.