

CaMKII Rabbit mAb



ABZ4888


Product Overview

Name	CaMKII Rabbit mAb
Accession(Primary)	Q13554
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 45,60 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200


Target information (Q13554)

Gene Name	CAMK2B
About Gene Symbol	





Immunohistochemistry analysis of paraffin-embedded mouse colon, using CaMKII antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescent analysis of CaMKII in PC-12 using CaMKII antibody.