

ILK Rabbit mAb


ABZ4894

Product Overview


Name	ILK Rabbit mAb
Accession(Primary)	Q13418
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q13418)


Gene Name	ILK
About Gene Symbol	ILK (Integrin Linked Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ILK in K562 lysates using Integrin Linked ILK antibody.



Immunohistochemistry analysis of paraffin-embedded Human stomach using ILK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Integrin Linked ILK in 293 using ILK antibody.

