

## ILK Rabbit mAb


ABZ4895

### Product Overview

|                             |                                                        |
|-----------------------------|--------------------------------------------------------|
| <b>Name</b>                 | ILK Rabbit mAb                                         |
| <b>Accession(Primary)</b>   | Q13418                                                 |
| <b>Host Species</b>         | Rabbit                                                 |
| <b>Clonality</b>            | Monoclonal                                             |
| <b>Isotype</b>              | IgG                                                    |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                      |
| <b>Calculated MW</b>        | Calculated MW: 51 kDa; Observed MW: 51 kDa             |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; ICC/IF-1:200-1:500;<br>FC-1:50-1:100 |

### Target information (Q13418)

|                          |                                                                                                                                                                                                                                                                                                        |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | ILK                                                                                                                                                                                                                                                                                                    |
| <b>About Gene Symbol</b> | ILK (Integrin Linked Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Integrin linked ILK in Jurkat, Raw264.7, PC-12, K562 lysates using ILK antibody.