

ASH2L Rabbit mAb


ABZ4896

Product Overview


Name	ASH2L Rabbit mAb
Accession(Primary)	Q9UBL3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:200-1:1000

Target information (Q9UBL3)

Gene Name	ASH2L
About Gene Symbol	



Western blot analysis of ASH2L in K562 lysates using ASH2L antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using ASH2L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.