

PAK2 Rabbit mAb


ABZ4898

Product Overview

Name	PAK2 Rabbit mAb
Accession(Primary)	Q13177
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20; IP-1:20-1:50

Target information (Q13177)

Gene Name	PAK2
About Gene Symbol	PAK2 (P21 (RAC1) Activated Kinase 2) is a serine/threonine protein kinase involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PAK2 in HeLa lysates
using PAK2 antibody.