

TDP43 Rabbit mAb


ABZ4899

Product Overview

Name	TDP43 Rabbit mAb
Accession(Primary)	Q13148
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:500-1:1000; FC-1:200-1:1000; IP-1:20-1:50

Target information (Q13148)

Gene Name	TARDBP
About Gene Symbol	TARDBP (TAR DNA Binding Protein) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TDP43 in K562 lysates
using TDP 43 antibody.