

FOXO1 Rabbit mAb


ABZ4900

Product Overview


Name	FOXO1 Rabbit mAb
Accession(Primary)	Q12778
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (Q12778)


Gene Name	FOXO1
About Gene Symbol	FOXO1 (Forkhead Box O1) is a transcription factor involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of FoxO1a in Jurkat, THP-1, SKOV-3, MCF-7 and HeLa lysates using FOXO1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain carcinoma using FoxO1a antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of FOXO1 in HeLa using FoxO1a antibody.