

RNA Helicase A Rabbit mAb


ABZ4902

Product Overview

Name	RNA Helicase A Rabbit mAb
Accession(Primary)	Q08211
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 141 kDa; Observed MW: 141 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:20

Target information (Q08211)

Gene Name	DHX9
About Gene Symbol	



Western blot analysis of RNA Helicase A in Jurkat lysates using RNA Helicase A antibody.