

SET Rabbit mAb


ABZ4917

Product Overview

Name	SET Rabbit mAb
Accession(Primary)	Q01105
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:500-1:1000; FC-1:10-1:100; IP-1:20-1:50

Target information (Q01105)

Gene Name	SET
About Gene Symbol	



Western blot analysis of SET in HepG2 lysates
using SET antibody.