

CDK6 Rabbit mAb


ABZ4920

Product Overview


Name	CDK6 Rabbit mAb
Accession(Primary)	Q00534
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:100-1:1000

Target information (Q00534)


Gene Name	CDK6
About Gene Symbol	CDK6 (Cyclin Dependent Kinase 6) is a serine/threonine protein kinase involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CDK6 in HeLa lysates using CDK6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using CDK6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of CDK6 in K562 using CDK6 antibody.

