

Cytochrome C Rabbit mAb


ABZ4921

Product Overview


Name	Cytochrome C Rabbit mAb
Accession(Primary)	P99999
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P99999)


Gene Name	CYCS
About Gene Symbol	



Western blot analysis of Cytochrome C in Human kidney lysates using Cytochrome C antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using Cytochrome C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Cytochrome C in Hela using Cytochrome C antibody.