

FOXO4 Rabbit mAb


ABZ4922

Product Overview


Name	FOXO4 Rabbit mAb
Accession(Primary)	P98177
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; FC-1:50-1:100

Target information (P98177)

Gene Name	FOXO4
About Gene Symbol	FOXO4 (Forkhead Box O4) is a transcription factor involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of FoxO4 in Human heart lysates using FOXO4 antibody.



Immunofluorescence analysis of FOXO4 in Hela cells using FoxO4 antibody.