

YB1 Rabbit mAb


ABZ4923

Product Overview


Name	YB1 Rabbit mAb
Accession(Primary)	P67809
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:200; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P67809)


Gene Name	YBX1
About Gene Symbol	



Western blot analysis of YB1 in HeLa lysates using YB1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using YB1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of YB1 in MCF-7 using YB1 antibody.

