

SUMO1 Rabbit mAb


ABZ4925

Product Overview


Name	SUMO1 Rabbit mAb
Accession(Primary)	P63165
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 80,12 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:20-1:100; IP-1:20

Target information (P63165)

Gene Name	SUMO1
About Gene Symbol	SUMO1 (Small Ubiquitin Like Modifier 1) is a metabolic enzyme involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SUMO1 in HeLa lysates using SUMO1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Sumo 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.