

Ubiquitin 60S Ribosomal Protein L40 Rabbit mAb


ABZ4927

Product Overview


Name	Ubiquitin 60S Ribosomal Protein L40 Rabbit mAb
Accession(Primary)	P62987
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 10 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:200-1:1000

Target information (P62987)


Gene Name	UBA52
About Gene Symbol	



Western blot analysis of UBA52 in C6 lysates using UBA52 antibody.



Western blot analysis of UBA52 in HeLa, SKOV-3 lysates using UBA52 antibody.



Immunocytochemistry analysis of UBA52(green)
in HeLa using UBA52 antibody, and DAPI(blue)