

Lysozyme Rabbit mAb


ABZ4929

Product Overview


Name	Lysozyme Rabbit mAb
Accession(Primary)	P61626
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:5000; ICC/IF-1:100-1:200

Target information (P61626)

Gene Name	LYZ
About Gene Symbol	



Western blot analysis of Lysozyme in K562 lysates; using Lysozyme antibody.



Immunohistochemistry analysis of paraffin-embedded mouse kidney using Lysozyme antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.