

RhoA/B/C Rabbit mAb


ABZ4931

Product Overview


Name	RhoA/B/C Rabbit mAb
Accession(Primary)	P61586
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:200-1:500; FC-1:200-1:500

Target information (P61586)


Gene Name	RHOA/RHOB/RHOC
About Gene Symbol	



Western blot analysis of RhoA/B/C in Jurkat lysates using RhoA/B/C antibody.



Immunocytochemistry analysis of RhoA/B/C(green) in HeLa using Rho A/B/C antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using RhoA/B/C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.