

Islet 1 Rabbit mAb


ABZ4933

Product Overview


Name	Islet 1 Rabbit mAb
Accession(Primary)	P61371
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:2000; ICC/IF-1:50-1:100; FC-1:20; IP-1:20

Target information (P61371)

Gene Name	ISL1
About Gene Symbol	



Immunocytochemistry analysis of Islet 1 (green) in SH-SY5Y using Islet 1 antibody, and DAPI (blue).



Western blot analysis of Islet 1 in Jurkat, K562 lysates using Islet 1 antibody.