


PTEN Rabbit mAb

ABZ4934


Product Overview

Name	PTEN Rabbit mAb
Accession(Primary)	P60484
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P60484)

Gene Name PTEN**About Gene Symbol**

Western blot analysis of PTEN in Raw264.7
lysates using PTEN antibody.



Western blot analysis of PTEN in MCF-7 lysates
using PTEN antibody.