

FOXA1 Rabbit mAb


ABZ4935

Product Overview

Name	FOXA1 Rabbit mAb
Accession(Primary)	P55317
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:20; ICC/IF-1:20; FC-1:20

Target information (P55317)

Gene Name	FOXA1
About Gene Symbol	FOXA1 (Forkhead Box A1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of FOXA1 in HepG2 lysates using FOXA1 antibody.