

VCP Rabbit mAb


ABZ4936

Product Overview

Name	VCP Rabbit mAb
Accession(Primary)	P55072
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 89 kDa; Observed MW: 89 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:10-1:100

Target information (P55072)

Gene Name	VCP
About Gene Symbol	



Western blot analysis of VCP in HeLa lysates using VCP antibody.