

NUP98 Rabbit mAb

ABZ4939


Product Overview

Name	NUP98 Rabbit mAb
Accession(Primary)	P52948
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 198 kDa; Observed MW: 188 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (P52948)

Gene Name NUP98

About Gene Symbol



Western blot analysis of NUP98 in HeLa lysates using NUP98 antibody.



Immunohistochemistry analysis of paraffin-embedded Human testis using NUP98 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.