

USP11 Rabbit mAb


ABZ4940

Product Overview

Name	USP11 Rabbit mAb
Accession(Primary)	P51784
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 110 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:200-1:1000; IP-1:10-1:100

Target information (P51784)

Gene Name	USP11
About Gene Symbol	USP11 (Ubiquitin Specific Peptidase 11) is a deubiquitinase involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and



Western blot analysis of USP11 in (1) HEK293 lysates; (2) mouse testis lysates using USP11 antibody.