


CCR3 Rabbit mAb


ABZ4941

Product Overview

Name	CCR3 Rabbit mAb
Accession(Primary)	P51677
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (P51677)

Gene Name CCR3**About Gene Symbol**



Immunohistochemistry analysis of paraffin-embedded Human epidermal of skin using CCR3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.