

DUSP3 Rabbit mAb


ABZ4942

Product Overview


Name	DUSP3 Rabbit mAb
Accession(Primary)	P51452
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20; FC-1:50-1:100; ICC/IF-1:50-1:100

Target information (P51452)


Gene Name	DUSP3
About Gene Symbol	



Western blot analysis of DUSP3/VHR in HeLa, 3T3, C6, Jurkat lysates using DUSP3 antibody.



Immunocytochemistry analysis of DUSP3(green) in HeLa using DUSP3 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using DUSP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.