

Rab7 Rabbit mAb


ABZ4943

Product Overview


Name	Rab7 Rabbit mAb
Accession(Primary)	P51149
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:500; FC-1:10-1:100

Target information (P51149)


Gene Name	RAB7A
About Gene Symbol	RAB7A (RAB7A, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of RAB7 in (1) A375 lysates;
(2) 3T3 lysates using Rab7 antibody.



Immunohistochemistry analysis of paraffin-
embedded Human colon using RAB7 antibody. High-
pressure and temperature Sodium Citrate pH 6.0
was used for antigen retrieval.



Immunofluorescence analysis of Rab7 in HeLa
using RAB7 antibody.

