

Mre11 Rabbit mAb

ABZ4945


Product Overview

Name	Mre11 Rabbit mAb
Accession(Primary)	P49959
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P49959)

Gene Name MRE11

About Gene Symbol



Western blot analysis of Mre11 in HeLa, K562, Jurkat lysates using Mre11 antibody.

